

## TECHNICAL SPECIFICATIONS

Energy Source	Bipolar balanced output (proprietary wave form)
Output Impedance	10 kilo Ohm
Outputs	2 independent outputs: 1 Face treatment hand piece 1 Body Treatment hand piece
Power Output Adjustment	Linear, continues adjustment from 0 to maximum
Electrical	220 V, 50 c/s or 110 V, 60 c/s-50 V/A
Body Probe Diameter	60mm
Facial Probe Diameter	28mm
Dimensions	35cm (length) x 27cm (width) x 24cm (height)
Weight	5,0 kg /11 lbs

## ABOUT VIORA

Viora is a market leader in the development of technologically advanced medical aesthetic systems for physicians and aesthetic clinicians worldwide. Dedicated to the research and development of safe and effective aesthetic solutions, Viora provides the highest standard of professional grade solutions while keeping your business in mind.

Viora offers a wide range of applications including body contouring, cellulite reduction, skin tightening procedures as well as phototherapy and anti-aging treatments. All of the products are designed by a team of veteran medical device developers in close cooperation with the company's advisory board including medical aesthetics practitioners, plastic surgeons, and dermatologists. All of Viora's products are manufactured under ISO certification No. 13485.

## CORPORATE HEADQUARTERS

30 Montgomery St. Suite 660  
Jersey city, NJ 07302  
U.S.A.  
T: +1.201.332.4100  
F: +1.201.332.4200  
Toll free: +1.877.384.8032  
Email: [info@vioramed.com](mailto:info@vioramed.com)

[WWW.VIORAMED.COM](http://WWW.VIORAMED.COM)



  
**VIORA**  
Revealing your potential

**Infusion™**  
Revealing healthy, youthful skin



The new **wave**  
of **needle-free** mesotherapy

  
**VIORA**  
Revealing your potential

## INFUSION™ ELECTRO-MESOTHERAPY THE PAIN-FREE, RISK-FREE ALTERNATIVE


Infusion™ Electro-mesotherapy is the safe, no needle alternative to traditional mesotherapy. It offers all the benefits of mesotherapy with none of the risks. Infusion™ revitalizes the skin with enriching solutions that leave it feeling supple, refreshed, and radiant. Based on the well established science of electroporation, Infusion™ maximizes cell permeability and absorption and enhances the overall appearance of the face and body.

Designed for aesthetic professionals, Infusion™ offers the perfect, non-invasive solution to a wide variety of cosmetic conditions and significantly expands treatment options with such aesthetic mainstays as:

- Anti aging and wrinkle reduction
- Collagen regeneration
- Skin rehydration
- Skin firming and toning
- Body sculpting and cellulite reduction
- Stretch mark reduction

## REFINING THE SCIENCE OF NEEDLE-FREE MESOTHERAPY

The Infusion™ Ionwave™ increases skin permeability to allow the safe passage of natural healing compounds. Precisely crafted pulse waves temporarily open micro channels through the cell membrane to bypass the stratum corneum layer without puncturing the skin.



### Ionwave™ Effects:

1. When the Ionwave™ properties come in contact with the cell they begin the structural rearrangement of the lipid bilayer to form micro pores in the cell membrane.



2. As the gatekeepers to the cell, micro channels enable or restrict the passage of natural nutritional compounds and extracts. Micro channels are short-lived and reseal soon after compounds are delivered.



## INFUSION™ IONIC SOLUTIONS STERILE, HIGH-QUALITY, HYPOALLERGENIC AND NON-IRRITATING



All natural, Viora solutions are created from the purest, high quality ingredients to maximize cell absorption and ensure safety and effectiveness. Scientifically engineered to achieve the desired effect, the interplay between the Ionic solution and the Infusion™ device determines the penetration depth and ensures sufficient permeability. Each solution is a blend of special serums, proteins and healing agents and is precisely formulated for the specific application. The use of these unique Infusion™ solutions in conjunction with the Infusion™ system is highly recommended to attain optimal results.

## INFUSION™- ENHANCE YOUR SERVICE PORTFOLIO

The Infusion™ system brings the highest quality aesthetic device to your clinic. It can be used alone or together with other treatments and is the perfect complement to any standard aesthetic procedure such as peeling and exfoliation. The optimal addition for the aesthetic practice, Infusion™ offers a comfortable and convenient solution for all your patient's skincare needs.

## INFUSION™ ADVANTAGES:

- Ionwave™: State-of-the-art technology designed to deliver optimal permeability without risk of injury or risk of infection.
- Immediately visible results: Maximum absorption of custom ionic solutions for results patients see and feel.
- No needles: Patients enjoy safe, painless treatments without the bruising, edema or infections.
- Intelligent, ergonomic design: Easy-to-use and functional device allows practitioners to comfortably perform treatments.
- Compact and lightweight unit: Easily portable for full treatment flexibility.



**Infusion™**  
Revealing healthy, youthful skin